

A digital microfluidics-based electrochemical impedance spectroscopy for cell-based immunoassay detection in a dynamic mode

Yuqian Zhang, and Yuguang Liu

Department of Surgery, Division of Surgical Research, Mayo Clinic, USA
Microbiome Program, Center for Individualized Medicine, Mayo Clinic, USA and
Department of Immunology, Mayo Clinic, USA

Presenter: Yuqian Zhang

MicroTAS 2022


Oct. 25th, 2022

Human peripheral blood mononuclear cells (PBMCs)

- Key drivers of the immune responses to pathogens
- Heavily involved in physiological homeostasis regulatory mechanisms
- immune-mediated indicators of therapeutic responsiveness

Advantages as biomarker:


- Minimally invasive
- Provide a more comprehensive overview of immune status



<https://www.the-scientist.com/>


Dynamic quantification of PBMC abundance

Standard method: **flow cytometry**



- Sensitive and reliable
- Relies on highly specialized and bulky equipment
- Not realistic to monitor PBMC abundance on a regular basis

Our method : **digital microfluidics (DMF)-based electrochemical impedance spectroscopy (EIS)**



DMF device

Dropbot system

- Palm-sized device, automated and parallel operation
- Cost-effective, minimal human intervention
- Ideally for point-of-care (POC) testing

Liu, Y., et al. *ACS omega* 6.39 (2021): 25642-25651.


Digital microfluidics (DMF)

Conventional lab-on-chip:

- Permanently etched microchannels
- External modules: pumps and microvalves

Digital microfluidic lab-on-chip:

- Manipulation of liquids as discrete droplets
- Electrowetting-on-dielectric (EWOD):



- Automated and programmable, reconfigured on-demand,



Liu, Y., et al. *ACS omega* 6.39 (2021): 25642-25651.


Electrochemical impedance spectroscopy (EIS)

Electronic-based detection




- Easy-to-miniaturize, rapid and real-time detection

Impedance (EIS) biosensor

- Utilizes a small amplitude, AC signal to probe the impedance characteristics
- Non-faradaic: change in the **double layer capacitance (C_{dl})** element
- (no direct current needed): **minimum destructive on electrodes, label-free detection**





Layout of the integrated DMF device



Overview of DMF device:

- **bottom plate:** actuation electrode array
- **top plate:** IDE sensing electrode


Equivalent Circuit and Modeling of IDEs



$$\frac{1}{|Z|} = \frac{1}{|Z_1|} + \frac{1}{|Z_2|}$$

$$|Z_1| = \sqrt{R_{sol}^2 + \frac{1}{(2\pi f C_{dl})^2}}$$

$$|Z_2| = \frac{1}{2\pi f C_{de}}$$




Two parallel branches:

Z_1 : capacitive and resistive responses of droplet solution



Z_2 : solution dielectric capacitance

Surface conjugation

Formation of thiol-carboxyl *self-assembled monolayer (SAM)*



PBMC hybridization:






Impedance change (%)

$$= \frac{Z_{cell} - Z_{base}}{Z_{base}} * 100$$



- 1) Blank control: 0.19%
- 2) non-specific binding: 2.2%
- 3) Target binding: 17.6%

Effect of IDE geometry



Over 3-fold sensitivity improvement

Dynamic incubation mode: *forward-backward-backward-forward*



Comparison of impedance increment (dynamic vs stationary modes)

PBMC Concentration (#/mL)	10^4	5×10^4	10^5	10^6
Impedance Increment (%)	242.7%	64.4%	26.9%	12.9%

- A novel DMF platform integrated with an EIS-based biosensor for the detection of PBMCs
- **Low sample volume (4 μ L) for rapid detection (20 min)**
- *Quantitative* detection of PBMC abundance **dynamic incubation modes showed 2.4-fold enhanced detection signal** detect as low as **10^4 PBMCs/mL**, approximately *two orders of magnitude less* than the biologically relevant range

Acknowledgements



Alexander Revzin, Ph.D.
Consultant
Department of Physiology
& Biomedical Engineering
Mayo Clinic, USA



Seth Hara, Ph.D.
Senior engineer
Division of Engineering
Mayo Clinic, USA



**Marina R. Walther-
Antonio, Ph.D.**
Associate Consultant I
Department of Surgery
Mayo Clinic, USA



Diane Eaker, M.S.
Senior engineer
Electronics Development
Unit of the Division of
Engineering

Funding resources: Ivan Bowen Family
Foundation, Mayo Clinic



- Department of Surgery,
- Microbiome Program,
Center for Individualized
Medicine

Thank you

Any questions?